



Entity
Medical



Read Me

ENTITY MEDICAL

Import , Supply Of Medical
Devices And Products



May Towers, Nasr City, Cairo 010 5583 4363 support@entitymedicalegypt.com



1688 AIM STRYKER CAMERA PLATFORM

Stryker's 1688 AIM 4K Platform is the first system to 4K fluorescence! The entire tower is managed by the Connected OR Hub, giving you seamless connectivity to images, data, and settings you need. Now, surgical teams can see in greater detail than ever with superior image clarity, more accurate color reproduction, and enhanced light quality. With the added benefit of 4K fluorescence

Brilliant Visualization

Designed from the ground up for seamless standardization and ease of use across multiple specialties, the 1688 AIM 4K Platform redefines surgical technology and patient care.

1688 AIM 4K Camera

Brilliant 4K resolution offers our surgeons crystal-clear images for improved visualization during surgery.

L11 light source and auto-light technology

The proprietary auto-light technology is designed to automatically adjust based on anatomy to provide consistent lighting, enabling our surgeons to focus on the procedure at hand and nothing else.

4K surgical display

A 32" monitor optimizes color in advanced imaging modes to allow surgeons to seamlessly transition between modes.

PneumoClear

This 50 liter insufflator is the first all-inclusive integrated heating, humidification, and smoke evacuation platform on the market. It provides consistently clear images through an inflow and outflow management system, designed for stability while actively evacuating smoke.

SPY-PHI

The 1688 AIM 4K Platform also supports SPY-PHI, a portable fluorescence imaging device that enables perfusion assessment in vessels and related tissue perfusion during plastic, microsurgical,

microsurgical, reconstructive and gastrointestinal procedures. The ability to visualize blood flow during open surgery may also help our surgeons make critical decisions that can reduce the risk of complications related to poorly perfused tissue.¹

1588 AIM STRYKER CAMERA PLATFORM

THE 1588 AIM IS OUR NEXT GENERATION VISUALIZATION PLATFORM INCORPORATING MULTIPLE ADVANCED IMAGING MODALITIES THAT ENABLE SURGEONS TO SEE AND DO MORE.

32" 4K SURGICAL DISPLAY

4X THE RESOLUTION OF A 1080P IMAGE WITH 4K PIXEL DENSITY, OFFERING A BRIGHT, COLORFUL VIEWING EXPERIENCE.

SDC3 DIGITAL CAPTURE DEVICE

STERILE FIELD DEVICE CONTROL, PHYSICIAN SPECIFIC PROFILES, WIFI CONNECTIVITY, AND CRITICAL DEVICE OVERLAY. EXPORTING OPTIONS INCLUDE: DVD, USB.

PNEUMOCLEAR INSUFFLATOR

INTEGRATED SMOKE EVACUATION, HEATING, AND HUMIDIFICATION FEATURES FOR A CONSISTENTLY CLEAR IMAGE AND IMPROVED PATIENT RECOVERY.

CLARITY IMAGE ENHANCEMENT SYSTEM

IMPROVES VISUALIZATION BY MITIGATING THE IMPACT SURGICAL PLUME, FOG, CONDENSATION, AND MURKY FLUID CAN HAVE ON THE IMAGE CLARITY.

1588 AIM 3-CHIP CAMERA CONTROL UNIT

IMPROVED IMAGE CLARITY, ADVANCED IMAGING AND STANDARDIZATION IN NINE SURGICAL SPECIALTIES.

L10 AIM LIGHT SOURCE

TRANSMITS WHITE LIGHT, NEAR-INFRARED, AND INFRARED WAVELENGTHS FOR MULTIPLE ADVANCED IMAGING CAPABILITIES.



1488 - 1288 - 1188

CAMERA STRYKER

• 1488

- STANDARD DESIGN WITH REMOVABLE C-MOUNT
COUPLER - INTEGRATED CAMERA AND COUPLER DESIGN -
PENDULUM/UROLOGY CAMERA DESIGN - STANDARD DESIGN WITH
OPTICAL ZOOM COUPLER



• 1288

DIGITAL/ANALOG: TWO DIGITAL VIDEO INTERFACE (HDTV)
FORMAT ENHANCEMENT: 16 LEVELS (SWITCHABLE) IMAGING
SYSTEM: 1/3" PROGRESSIVE SCAN CCDS, HIGH DEFINITION



• 1188


- STANDARD CABLE DESIGN FOR MULTI-SPECIALTY USE -
TEXTURED FINISH FOR SECURE HANDLING - CAMERA BUTTON
OFFERS 6 FUNCTIONS FOR IMPROVED SURGEON CONTROL -



L71 LED LIGHT SOURCE
with Auto lightTechnology



FIBEROPTIC LIGHT CABLES



VISION PRO LED DISPLAY



L9000 LED LIGHT SOURCE

CORE Console

Powered Instrument Driver

The Stryker CORE Console Energized Implement Operator seamlessly empowers surgical procedures with its sophisticated built-in fluid propulsion mechanism. This cutting-edge device incorporates a streamlined touchscreen interface for facile user interaction, allowing for tailoring of preferences and seamless control of flow rates. Additionally, it features an integrated support system that enhances procedural efficacy, marking it as an advanced solution for surgical interventions.



CORE 2 Console

The CORE 2 Console is Stryker's next gen console that enhances surgical workflows with customizable user profiles, smart programming, and integrated irrigation. The console supports a wide range of surgical applications and with its intuitive interface, safety features, and ability to power multiple devices makes it a reliable and cost-effective choice for healthcare professionals.



Formula Shaver Handpieces

Resection


Stryker's CrossBlade and Formula shaver lines offer varying working lengths, window sizes and window styles to accommodate many arthroscopic procedures.



L71 LED LIGHT SOURCE
with Auto lightTechnology




FIBEROPTIC LIGHT CABLES



VISION PRO LED DISPLAY



L9000 LED LIGHT SOURCE



4K 32" OLED MONITOR



INSUFFLATOR 40L



PNEUMOSURE XL INSUFFLATOR



PNEUMOCLEAR INSUFFLATOR

1488 - 1288 - 1188

CAMERA STRYKER

• 1488

- STANDARD DESIGN WITH REMOVABLE C-MOUNT
COUPLER - INTEGRATED CAMERA AND COUPLER DESIGN -
PENDULUM/UROLOGY CAMERA DESIGN - STANDARD DESIGN WITH
OPTICAL ZOOM COUPLER



• 1288

DIGITAL/ANALOG: TWO DIGITAL VIDEO INTERFACE (HDTV)
FORMAT ENHANCEMENT: 16 LEVELS (SWITCHABLE) IMAGING
SYSTEM: 1/3" PROGRESSIVE SCAN CCDS, HIGH DEFINITION




• 1188

- STANDARD CABLE DESIGN FOR MULTI-SPECIALTY USE -
TEXTURED FINISH FOR SECURE HANDLING - CAMERA BUTTON
OFFERS 6 FUNCTIONS FOR IMPROVED SURGEON CONTROL -



L71 LED LIGHT SOURCE
with Auto lightTechnology



VISION PRO LED DISPLAY



FIBEROPTIC LIGHT CABLES



L9000 LED LIGHT SOURCE

CORE Console

Powered Instrument Driver

The Stryker CORE Console Energized Implement Operator seamlessly empowers surgical procedures with its sophisticated built-in fluid propulsion mechanism. This cutting-edge device incorporates a streamlined touchscreen interface for facile user interaction, allowing for tailoring of preferences and seamless control of flow rates. Additionally, it features an integrated support system that enhances procedural efficacy, marking it as an advanced solution for surgical interventions.



CORE 2 Console

The CORE 2 Console is Stryker's next gen console that enhances surgical workflows with customizable user profiles, smart programming, and integrated irrigation. The console supports a wide range of surgical applications and with its intuitive interface, safety features, and ability to power multiple devices makes it a reliable and cost-effective choice for healthcare professionals.




Formula Shaver Handpieces

Resection

Stryker's CrossBlade and Formula shaver lines offer varying working lengths, window sizes and window styles to accommodate many arthroscopic procedures.



L71 LED LIGHT SOURCE
with Auto lightTechnology




VISION PRO LED DISPLAY



FIBEROPTIC LIGHT CABLES



L9000 LED LIGHT SOURCE



4K 32" OLED MONITOR



INSUFFLATOR 40L



PNEUMOSURE XL INSUFFLATOR



PNEUMOCLEAR INSUFFLATOR

IMAGE 1 S STORZ

26" FULL HD MONITOR . IMAGE1 S CONNECT .

USB ADAPTOR, FOR ACC 1 AND ACC 2 .

IMAGE1 S H3-LINK .

IMAGE1 S H3-Z FI THREE-CHIP

FULL HD CAMERA HEAD . COLD LIGHT

FOUNTAIN XENON 300 SCB .

ENDOFLATOR® 50 SCB .

HAMOU®/ENDOMAT® SCB . AUTOCON® 200 .

TWO-PEDAL FOOTSWITCH USB .

NEUTRAL ELECTRODE . NEUTRAL ELECTRODE

CONNECTING CABLE .

TWO-PEDAL FOOTSWITCH .

AIDA™ DOCUMENTATION SYSTEM .

EQUIPMENT CART, WIDE .



**CAMERA karl storz
image one hub**



**CAMERA karl storz
image one spies**



**CAMERA karl storz
image one hub
with usb**



IMAGE 1 SPIES STORZ

CREATES INNOVATION!

KARL STORZ IMAGE1 SPIES ONCE AGAIN SETS NEW STANDARDS

IN IMAGING FOR MINIMALLY INVASIVE SURGERY. THE IMAGE1 SPIES MODULAR DESIGN ALLOWS CUSTOMIZATION

TO INDIVIDUAL NEEDS SO THAT THE SYSTEM IS OPTIMALLY ADAPTED TO SUIT USER REQUIREMENTS. THE SPIES

TECHNOLOGY PROVIDES VIDEO ENDOSCOPES WITH DIVERSE VISUALIZATION POSSIBILITIES FOR ENDOSCOPIC

IMAGING. THE MODULAR CONCEPT OF IMAGE1 SPIES INTEGRATES TODAY'S VIDEO ENDOSCOPES FROM KARL

STORZ INTO THE NEW CAMERA PLATFORM VIA IMAGE1 X-LINK™



**TH 100
H3 - Z S**



**TH 102
H3 - Z FI S**



**TH 103
H3 - P S**



**TH 104
H3 - ZA S**



**TH 120
4K - UHD**



**TH 110
HX**

The Valleylab™ LS10

Valleylab LigaSure

DIMENSIONA AND WEIGHT »

WIDTH 300 MM (11.81 INCHES)

» DEPTH 377 MM (14.84 INCHES)

» HEIGHT 105 MM (4.13 INCHES)

» WEIGHT 5 KG (11 LBS)

AUDIO VOLUME

THE STATED AUDIO LEVELS ARE AT A DISTANCE OF ONE METER.

ALERT TONES MEET THE REQUIREMENTS OF IEC60601-2-2.



The Valleylab™ LS8

Valleylab LigaSure

DIMENSIONA AND WEIGHT

» WIDTH: 38.6 CM (15.2 IN)

» DEPTH: 40.6 CM (16.0 IN)

» HEIGHT: 12.7 CM (5.0 IN)

» WEIGHT: 5.9 KG (13 LBS) TYPICAL

AUDIO VOLUME

THE STATED AUDIO LEVEL IS FOR THE ACTIVATION TONE AND ALARM TONE AT A DISTANCE OF ONE METER. ALARM TONES MEET THE REQUIREMENTS OF IEC 60601-2-2.



SCOPE

HYSTEROSCOPY SET (INSPECTION)




HYSTEROSCOPY SET (WORK)



HYSTERO RESECTOSCOPE



VAPORIZAION RESECTOSCOPE



CYSTOSCOPE (LOCKABLE)



OLYMPUS CLV-190

One-touch access to HDTV and NBI

MAIN FEATURES

AUTOMATIC LIGHT ADJUSTMENT TO ACHIEVE THE IDEAL ILLUMINATION FOR OBSERVATION WITH EACH SCOPE



NBI (NARROW BAND IMAGING) IN EVIS EXERA III 190 SERIES SCOPES PROVIDES TWICE THE VIEWABLE DISTANCE AND A SIGNIFICANTLY BRIGHTER IMAGE, THANKS TO AN IMPROVED LAMP DESIGN AND SIGNAL PROCESSING.

THE NEWLY DESIGNED, WATERPROOF ONE-TOUCH CONNECTOR ALLOWS A ONE-STEP CONNECTION TO THE LIGHT SOURCE AND DOES NOT REQUIRE A SCOPE CABLE.

CONSIDERABLE REDUCTION IN OPERATING NOISE, THANKS TO REDESIGNED FAN CONNECTION TO THE LIGHT SOURCE AND DOES NOT REQUIRE A SCOPE CABLE.

KARL STORZ

ENDO-SCOPE

SPECIAL FEATURES

- DISTAL END OF SHEATH AUTOMATICALLY SHAPED, WITH ROUNDED TIP .
- MINIMAL SHEATH DIAMETER .
- THE LARGE WORKING CHANNEL ALLOWS THE USE OF RIGID INSTRUMENTS AND PROBES UP TO 5 FR .
- MAXIMUM INFLOW DUE TO LARGE CENTRAL WORKING CHANNEL WITH 2 LATERAL ,RIGHT-AN- GLED IRRIGATION PORTS .

URETERO-RENSCOPE,AUTOCLAVABLE,
LENGTH 43 CM
URETERO-RENSCOPE,AUTOCLAVABLE,
LENGTH 34 CM

SPECIFICATIONS

- DISTAL TIP
- INSTRUMENT SHEATH
- WORKING CHANNEL
- TELESCOPE
- LENGTH
- EYEPIECE



CLV-180

EVIS EXERA II Xenon Light Source

AUTOMATICALLY ADJUSTS LIGHT INTENSITY TO ACHIEVE IDEAL ILLUMINATION FOR OBSERVATION OF THE GASTROINTESTINAL TRACT



HIGH-INTENSITY 300 W LIGHT SOURCE DELIVERS THE ILLUMINATION TO DRIVE HDTV AND A NEW NARROW BAND IMAGING FEATURE ENHANCES VISUALIZATION CAPABILITIES.

EQUIPPED WITH SPECIALLY COATED FILTERS FOR NBI (NARROW BAND IMAGING)

COMPATIBLE WITH:

SURGICAL ENDOSCOPES (HD ENDOSCOPE AND ENDOSCOPE VIDEOSCOPIES, FLEXIBLE VISERA VIDEOSCOPIES, AND 1CCD/3CCD CAMERA HEADS) AND EVIS 100/130/140/160/180 SERIES ENDOSCOPES AS WELL AS BRONCHOSCOPES

TELESCOPES AND LAPAROSCOPES



KARL STORZ Arthroscope Sheaths

for fast and secure telescope-sheath connection



KARL STORZ arthroscope sheaths are characterized by the following special features:

STABLE AND SECURE TELESCOPE-SHEATH CONNECTION

1. STANDARD COUPLING MECHANISM
2. SNAP-IN COUPLING MECHANISM
3. AUTOMATIC LOCK-IN COUPLING MECHANISM

• HOMOGENEOUS MECHANICAL STABILITY

• HOMOGENOUS TRANSITION FROM THE OBTURATOR TO THE DISTAL TIP ENSURES ATRAUMATIC ENTRY OF THE ARTHROSCOPE SHEATH INTO THE JOINT

• SHEATHS CAN BE FULLY DISMANTLED FOR EFFECTIVE CLEANING AND REPROCESSING

• A SPECIALLY DESIGNED GASKET AT TACHMENT ENABLES RAPID CLOSURE OF THE SHEATHS AND THEIR USE AS A WORKING CHANNEL FOR INTRODUCING INSTRUMENTS. THIS GASKET ATTACHMENT CAN BE USED WITH ALL KARL STORZ ARTHROSCOPE SHEATHS

SYSTEM 7

Heavy Duty Power Tools

PERFORMANCE

◆ IMPROVED POWER-TO-WEIGHT RATIO RESULTS IN A PERFECT BALANCE OF SPEED, WEIGHT AND TORQUE

◆ ADVANCED MICROPROCESSOR CONTINUALLY MONITORS BOTH THE HANDPIECE AND BATTERY FOR OPTIMAL PERFORMANCE

◆ ROTARY HANDPIECE OFFERS SUBSTANTIALLY INCREASED

TORQUE, PROVIDING OPTIMUM PERFORMANCE

◆ RE-ENGINEERED SAGITTAL SAW BLADE MOUNT OFFERS SIGNIFICANT NOISE REDUCTION

Efficiency

◆ TRIATHLON ATTACHMENT ELIMINATES THE NEED FOR EXTRA ATTACHMENTS ◆ UNIQUE ROTARY HANDPIECE DESIGN ALLOWS ALL ATTACHMENTS TO BOTH DRILL AND REAM

EASE OF USE

• IMPROVED ERGONOMICS:

• 29% REDUCTION IN OVERALL


SYSTEM WEIGHT

• IMPROVED SAGITTAL HEAD DESIGN OFFERS SIMPLE, SECURE BLADE LOADING

• SYSTEM DESIGNED AND VALIDATED FOR AUTOMATED WASHING
• ROTARY HANDPIECE TRIATHLON ATTACHMENT REDUCES DISTANCE TO THE SURGICAL SITE



ADDITIONAL SERVICES




BIPOLAR GRASPING FORCEPS

FLAT JAWS, SERRATED,
SIZE 2.4 MM, LENGTH 26 CM




MINIATURE STRAIGHT FORWARD TELESCOPE

SEMRIGID WITH REMOTE
EYEPIECE, AUTOCLAVABLE



HOPKINS FORWARD OBLIQUE TELESCOPE 0

ENLARGED VIEW,
DIAMETER 4MM, LENGTH
30CM, AUTOCLAVABLE



BIGATTI INTRAUTERINE SHAVER

PROCEDURE FOR MAJOR
HYSTEROSCOPIC OPERATIONS



HOPKINS FORWARD OBLIQUE TELESCOPE 30
ENLARGED VIEW,

DIAMETER 4MM, LENGTH
30CM, AUTOCLAVABLE

**Micro Alligator Forcep
Left, Cup Shaped**

SIZE, CUP SIZE, JAW SIZE
8 CM (3"), 1.0 X 0.9 MM, 4.0 MM
8 CM (3"), 0.6 X 0.5 MM, 3.5 MM



**Anderson-Adson
Self Retaining Retractor
4 x 4 Sharp Prongs**

SIZE
20 CM (8")
21.5 CM (8 ½")
22 CM (8 ¾")




ADDITIONAL SERVICES




SIZE
30.5 CM (12")
32 CM (12 ½")

**Beckmann-Adson
Self Retaining Retractor**



SIZE
13 CM (5")
14 CM (5 ½")
15.5 CM (6")

**Weitlaner-Wullstein
Self Retaining Retractor
3 x 3 Blunt Prongs**



SIZE
27.5 CM (10 ¾")
23.5 CM (9 ¼")
27.5 CM (10 ¾")
26 CM (10 ¼")
25.5 CM (10")
24.5 CM (9 ¾")

**Stille-Liston Bone
Cutting Forcep
Straight, Compound Action**

**Ferris-Smith Kerrison
Rongeur
Detachable Model, Up Cutting**

SIZE
18 CM (7"), 1 MM
18 CM (7"), 2 MM
18 CM (7"), 3 MM
18 CM (7"), 4 MM
18 CM (7"), 5 MM



**Lemnectomy Rongeur
Straight, Fenestrated and
Serrated Jaws**

SIZE, BITE SIZE, PROFILE
15.5 CM (6"), 4 X 14 MM, STRAIGHT
15.5 CM (6"), 4 X 14 MM, UP
15.5 CM (6"), 4 X 14 MM, DOWN
18 CM (7"), 5 MM

